

Integro Surgical Products

# Instruction Guide for Solution Stands



Read the instructions and safety information in this manual before operating this product.

## Table of Contents

Features.....	3
Care and Maintenance .....	4

## Unpacking and Inspection

InnerSpace solution stands are heavy. We recommended that two people assist in the unpacking.

- If the shipping container appears to be damaged in any way, contact the shipping company immediately. Save all damaged packaging materials.
- Carefully inspect the solution stand while unpacking. If any damaged or defective parts are seen, call InnerSpace Customer Service (888-435-2256).
- Solution stands are shipped fully assembled and ready for immediate use. No assembly is required.

# Features

- Fully welded stainless steel construction
- Single- and double-basin designs available
- Available with lower shelf
- Open frame design allows accessibility from any side
- 2" swivel casters

Product Number	Dimensions	Number of Basins	Basin Capacity	Lower Shelf
SMSSOL7QH	14" D x 14" W x 34" H	1	7 Quart	N
SMSSOL7QS	14" D x 14" W x 34" H	1	7 Quart	Y
SMDSOL7QH	28" D x 14" W x 34" H	2	7 Quart	N
SMDSOL7QS	28" D x 14" W x 34" H	2	7 Quart	Y



Stainless steel solution stands must be cleaned on a regular basis in order to prevent unnecessary damage to the stainless steel surface. When cleaning solution stands, make sure to use the proper cleaning agents and cleaning materials.

Do not use the following materials on the stainless steel surfaces of the solution stands:

- Abrasive pads
- Scrapers
- Steel wool
- Wire brushes

Do not use the following agents on the stainless steel surfaces of solution stands:

- Hard water
- Bleach or any products containing chlorine
- Hydrochloric acid
- Denatured alcohol

Always follow the manufacturer's directions for any chemicals used on the stainless steel surfaces of the solution stands.

Use the following cleaning materials and agents:

- Soft, clean lint-free cloth
- Distilled water
- Distilled water with mild detergent
- Cleaners approved for use on stainless steel

### **Cleaning a Stainless Steel Solution Stand**

- Using a damp, lint-free cloth and approved cleaner, wipe down all surfaces of the solution stand.
- Use a lint-free dry cloth to dry the solution stand or let air dry. If air drying, ensure solution stand is completely dry before storing. Store in a cool, dry place.

The solution stand can become damaged and unsafe if used improperly.

## **Solution Stands**

- Never overload the solution stand.
- Do not drop heavy loads on the solution stand.
- Do not use the solution stand at high speeds.

Using the solution stand in any of the above conditions can result in failure of the casters, wheels, or other parts of the equipment.

Periodically inspect the bottom of the solution stand and look for any of the following:

- Broken welds.
- Loose nuts or bolts. These should be tightened.
- Frame distortion that could be caused by an overloaded stand or by impacts with other objects while moving the stand.

## **Casters and Wheels**

- Make sure the casters are in good condition. If a swivel caster has excessive play due to wear, it may need to be replaced. If casters don't turn freely, corrosion or dirt may require their replacement.
- Wheels should be checked for visible tread wear. Flat spots due to accumulation of material, such as string or thread, can cause the wheel to bind. The axle bolt nut should be removed and the foreign material removed. If the internal components are intact, the caster can be reassembled.
- Rubber tires can be badly worn and cause erratic steering, bumping, load shifting, and floor damage.