

Integro Surgical Products

Instruction Guide for One-Step Step Stool



Read the instructions and safety information in this manual before operating this product.

Table of Contents

Features.....	3
Instructions.....	4
Care and Maintenance	6

Read these instructions before using. If step stools are used improperly, personal injury and damage to the step stools could result.

If components are damaged or missing, contact InnerSpace Customer Service (1-888-435-2256). Do not use substitute parts.

Safety Precautions

To prevent injury, please follow the precautions listed below before using the single step stool. Practicing these safety steps can greatly increase your odds of not being injured.

- Always make sure that the surface you place the step stool on is level. If all four feet are not stable, the risk of a fall-related injury is greater.
- Make sure the step stool is in the best position to perform the work. If you must reach out or bend excessively, stop and reposition the step stool. Additionally, stand near the middle of the stool as much as possible.
- Do not move the step stool while you are on it. Always descend the stool before repositioning it.
- Never attempt to increase the height of the step stool by placing it on top of boxes or other equipment.
- Inspect the step stool for damage regularly. Cracked steps can result in the total collapse of the step stool. Ensure that all four rubber floor tips are present.
- Do not use the step stool on wet, icy, or slippery surfaces.
- Children should not be left unsupervised around the step stool.


Standard

- Heavy-duty, stainless steel construction
- Tip resistant
- Corrugated, black rubber top surface
- Rubber floor tips
- Single-step stools are stackable and can connect together to form a large platform.

Optional

- Easy-roll transport cart with four swiveling casters (two locking). Available for single step stools (can hold up to five single step stools).

Product Number	Dimensions	Weight Capacity
SM121805STEP	12" D x 18" W x 6.25" H	650 lb



Connecting Single Step Stools to Create a Larger Platform

- Using the connecting hooks located on one long side and one short side of the single step stool, you can connect a number of these stools in a variety of configurations in order to form a larger working platform.
- All necessary safety precautions should still be observed when using the single step stools in any configuration.




Single step stool with connecting hooks/carry handles.



To hook step stools together, hook the lower lip of the long or short hookless side of one single step stool to the desired hook on a companion single step stool.

Instructions

Note the legs must be offset slightly when connecting along the long axis of both stools. Hooking the stools together with their legs adjacent will result in a wobbly and unsafe working platform.



Possible configurations of connected single step stools

Step stools must be cleaned on a regular basis in order to prevent unnecessary damage to the stainless steel surface. When cleaning the step stools, make sure to use the proper cleaning agents and cleaning materials.

Do not use the following materials on the stainless steel surfaces of the step stools:

- Abrasive pads
- Scrapers
- Steel wool
- Wire brushes

Do not use the following agents on the stainless steel surfaces of the step stools:

- Hard water
- Bleach or any products containing chlorine
- Hydrochloric acid
- Denatured alcohol

Always follow the manufacturer's directions for any chemicals used on the stainless steel surfaces of the step stools.

Use the following cleaning materials and agents:

- Soft, clean lint-free cloth
- Distilled water
- Distilled water with mild detergent
- Cleaners approved for use on stainless steel

Cleaning a Step Stool

- Using a damp, lint-free cloth, and approved cleaner, wipe all stainless steel interior and exterior surfaces of the step stool.
- Use any mild detergent to clean the corrugated rubber surfaces.
- Use a lint-free dry cloth to dry the step stool, or let it air dry. If air drying, ensure that the step stools are completely dry before storing, and store only in dry conditions.

Decals or printed labels:

- Use distilled water and a mild detergent applied with a lint-free cloth. Do not use any chemical cleaners.