

Ventaire™ Scope Drying and Tracking Cabinet



Support scope drying and storage processes, safety protocols, and reduce the chances (and time and expense) of reprocessing expired endoscopes with the Ventaire scope drying and tracking cabinet. The system meets the 2021 Multisociety Task Force Scope Storage and Drying Guidelines. The easy-to-operate system continuously filters air into internal scope channels to efficiently and thoroughly dry scopes and features timed drying and venting cycles. HEPA filtered air removes 99.9% retention rate of particles less than 0.1µm. The Ventaire system works with all major endoscope brands.

- Medical-grade plug and UL-approved electrical system
- HEPA filter removes 99.9% of particulates
- Energy-efficient, quiet air compressor flows air into manifold
- Manifold calibrates air circulation into active scope channels
- Filter controls air circulation into endoscopes
- Tubes connect to components in the Ventaire system to vent HEPA filtered air through endoscope channels
- Continuous venting of HEPA filtered air provides safe scope storage 24/7
- 10-minute default drying cycle
- Easy-to-operate 10" touchscreen
- Modify or add to scope inventory at any time
- Customizable report captures all scope information and status
- Monitor scope storage from touchscreen on cabinet or remotely via website (subscription required)
- Connect USB drive to USB port to generate report on scope activity within cabinet

AireCore Scope Cabinet (SEA1936SGS) with
8-port Ventaire Scope Drying and Tracking Unit (SVDUG8)


Ventaire Scope Drying and Tracking Cabinet




The 19.25"-deep Ventaire scope drying and tracking cabinet comes in 3 widths and 3 heights and with solid-front or tempered glass doors. It includes a Ventaire scope drying and tracking unit, fan ventilation unit, filter and tubes, and 10" touchscreen. Material choices include stainless steel, or nonporous AireCore with Glacier White finish, Brushed Aluminum finish, or specified laminate surface. The scope drying and tracking unit can be installed in an existing Evolve scope cabinet with a hinged door.

- HEPA fan unit provides positive-pressure air in cabinet
- Single and double adjustable scope managers adapt to different lengths of scopes. Single scope managers store 5, 8, or 9 scopes; double scope managers store 10, 16, or 18 scopes
- Medical-grade plug and UL-approved electrical system
- Manifold calibrates air circulation into active scope channels
- Energy-efficient, quiet air compressor flows air into manifold
- Adjustable cord and tube management keep endoscope tubes organized, secure, and elevated
- No lock, keyed, keyless, or proximity lock options
- Tubes connect to components in the Ventaire system to vent HEPA filtered air through endoscope channels
- Self-closing door
- Concealed hinge design
- Easy-to-operate 10" touchscreen
- Connect USB drive to USB port to generate report on scope activity within cabinet
- Padded rear wall protects scope optics
- Cabinets store scopes up to 78" in length
- Removable drip tray eases cleaning
- Air exits cabinet through vents
- Standard 4" h toe kick or optional 6" h toe kick


Ventaire Scope Drying and Tracking Cabinet




26.75\"/>




26.75\"/>



26.75\"/>



36\"/>



40\"/>

Dimensions (Nominal)

Exterior (interior), Scope hanging length:

26.75\"/>

26.75\"/>

26.75\"/>

36\"/>

40\"/>

Interior Configuration:

48\"/>

92\"/>

92\"/>

92\"/>

Options

- Solid-front and glass door on 92\"/>
- Door can be hinged on right or left side of 26.75\"/>
- Keyed, keypad, and proximity locks
- Sloped top
- 6\"/>
- Subscription to InnerSpace Smart Solutions website for remote monitoring of scopes

Specifications

Compressor air output:

- up to 45 LPM

Compressor decibels:

- 62 dB, with doors closed

Electrical requirements for drying unit:

- 120 VAC, 60Hz, 15-amp standard duplex receptacle
- Includes a hospital-grade, UL-approved, 125V, 10-amp, 10-foot-long power cord with Y-splitter to power both the drying and tracking unit and the fan unit

Electrical requirements for cabinet fan ventilation unit (SFVU):

- Voltage Rating: 115VAC
- Watts: 13.0W
- Amps: 0.19A

Ventaire Scope Drying and Tracking Cabinet

Replace these items every 6 months for optimal performance:

- Compressor tube/air filter replacement, package of 1: SVCTAF
- Manifold tube replacement, package of 5: SVMT5
- Manifold tube replacement, package of 8: SVMT8
- Manifold tube replacement, package of 9: SVMT9
- Manifold auxiliary water tube replacement, package of 5: SVWAT5
- Manifold auxiliary water tube replacement, package of 8: SVWAT8
- Manifold auxiliary water tube replacement, package of 9: SVWAT9
- HEPA filter replacement, package of 1: SHEPFIL, package of 3: SHEPFIL3


Single-use connector kits:

To find which connector kit is compatible with your endoscope, visit <https://www.innerspacehealthcare.com/Ventaire>

- Connector Kit 1, package of 100: SVCK1
Includes connector and biopsy port cover
- Connector Kit 2, package of 100: SVCK2
Includes connector, biopsy port cover, and auxiliary water connector
- Connector Kit 3, package of 100: SVCK3
Includes connector and luer plug
- Connector Kit 4, package of 100: SVCK4
Includes connector, water jet adaptor, and luer plug
- Connector Kit 5, package of 100: SVCK5
Includes connector, biopsy port cover, and luer coupler

To order, call an InnerSpace Customer Service representative at 888-435-2256 or visit <https://www.innerspacehealthcare.com/Ventaire>

Product Number Reference



Product Numbers

To order a Vентаire drying and tracking cabinet, first specify the cabinet size, door type, material, and configuration, and then specify a drying and tracking unit.

Scope Cabinet – 26.75"w x 19.25"d x 48"h

AireCore: SEA1927SGLW, SEA1927SGRW

AireCore with Brushed Aluminum: SEB1927SGLW, SEB1927SGRW

Stainless steel: SES1927SGLW, SES1927SGRW

Scope Cabinet – 26.75"w x 19.25"d x 60"h

AireCore: SEA1927SGLB, SEA1927SGRB

AireCore with Brushed Aluminum: SEB1927SGLB, SEB1927SGRB

Stainless steel: SES1927SGLB, SES1927SGRB

Scope Cabinet – 26.75"w x 19.25"d x 92"h

AireCore: SEA1927SGLS, SEA1927SGRS, SEA1927SSLS, SEA1927SSRS;
with double scope holder: SEA1927SGLD, SEA1927SGRD, SEA1927SSLD, SEA1927SSRD

AireCore with Brushed Aluminum: SEB1927SGLS, SEB1927SGRS, SEB1927SSLS, SEB1927SSRS;
with double scope holder: SEB1927SGLD, SEB1927SGRD, SEB1927SSLD, SEB1927SSRD

Stainless steel: SES1927SGLS, SES1927SGRS, SES1927SSLS, SES1927SSRS;
with double scope holder: SES1927SGLD, SES1927SGRD, SES1927SSLD, SES1927SSRD

Scope Cabinet – 36"w x 19.25"d x 92"h

AireCore: SEA1936SGS, SEA1936SSS;
with double scope holder: SEA1936SGD, SEA1936SSD

AireCore with Brushed Aluminum: SEB1936SGS, SEB1936SSS;
with double scope holder: SEB1936SGD, SEB1936SSD

Stainless steel: SES1936SGS, SES1936SSS; with double scope holder: SES1936SGD, SES1936SSD


Scope Cabinet – 40"w x 19.25"d x 92"h

AireCore: SEA1940SGS, SEA1940SSS;
with double scope holder: SEA1940SGD, SEA1940SSD

AireCore and Brushed Aluminum: SEB1940SGS, SEB1940SSS;
with double scope holder: SEB1940SGD, SEB1940SSD

Stainless steel: SES1940SGS, SES1940SSS; with double scope holder: SES1940SGD, SES1940SSD

Product Number Reference



Product Numbers

When ordering a drying and tracking unit for a stainless steel cabinet, add an S to the beginning of the product number.

Ventaire Scope Drying and Tracking Unit – For 26.75" w x 48" h left-hinged cabinet

5-port unit: SVDUGL5W (for glass door cabinet)

Ventaire Scope Drying and Tracking Unit – For 26.75" w x 48" h right-hinged cabinet

5-port unit: SVDUGR5W (for glass door cabinet)

Ventaire Scope Drying and Tracking Unit – For 26.75" w x 60" h left-hinged cabinet

5-port unit: SVDUGL5B (for glass door cabinet)

Ventaire Scope Drying and Tracking Unit – For 26.75" w x 60" h right-hinged cabinet

5-port unit: SVDUGR5B (for glass door cabinet)

Ventaire Scope Drying and Tracking Unit – For 26.75" w left-hinged cabinet

5-port unit: SVDUSL5 (for solid door cabinet), SVDUGL5 (for glass door cabinet)

10-port unit: SVDUSL10 (for solid door cabinet), SVDUGL10 (for glass door cabinet)

Ventaire Scope Drying and Tracking Unit – For 26.75" w right-hinged cabinet

5-port unit: SVDUSR5 (for solid door cabinet), SVDUGR5 (for glass door cabinet)

10-port unit: SVDUSR10 (for solid door cabinet), SVDUGR10 (for glass door cabinet)

Ventaire Scope Drying and Tracking Unit – For 36" w cabinet

8-port unit: SVDUS8 (for solid door cabinet), SVDUG8 (for glass door cabinet)

16-port unit: 16 port: SVDUS16 (for solid door cabinet), SVDUG16 (for glass door cabinet)

Ventaire Scope Drying and Tracking Unit – For 40" w cabinet

9-port unit: SVDUS9 (for solid door cabinet), SVDUG9 (for glass door cabinet)

18-port unit: SVDUS18 (for solid door cabinet), SVDUG18 (for glass door cabinet)

Ventaire Scope Drying and Tracking Cabinet

Materials and Finishes

Cabinet shell:

AireCore: Aluminum composite with Glacier White finish or specified Wilsonart, Formica, or Pionite standard laminate

AireCore with Brushed Aluminum finish: Aluminum composite with Brushed Aluminum finish

Stainless Steel: 18-gauge 304 stainless steel with 18% chrome and 8% nickel analysis; exposed surfaces

Edge banding:

AireCore with Brushed Aluminum finish: Brushed Aluminum acrylic

Solid-front door:

AireCore: Aluminum composite with Glacier White finish or specified Wilsonart, Formica, or Pionite standard laminate

AireCore with Brushed Aluminum finish: Aluminum composite with Brushed Aluminum finish

Stainless steel: 18-gauge 304 stainless steel with 18% chrome and 8% nickel analysis

Glass-front door:

AireCore: Aluminum composite with Glacier White finish or specified Wilsonart, Formica, or Pionite standard laminate with tempered glass and clear anodized aluminum trim

AireCore with Brushed Aluminum finish: Aluminum composite with Brushed Aluminum finish with tempered glass and clear anodized aluminum trim

Stainless steel: 18-gauge 304 stainless steel with 18% chrome and 8% nickel analysis and tempered glass

Hinges:

AireCore, and AireCore with Brushed Aluminum cabinet: Metal concealed European-style

Stainless steel cabinet: Stainless steel

Door pulls:

AireCore, and AireCore with Brushed Aluminum cabinet: Dull chrome-plated steel bow pull

Stainless steel cabinet: 5" stainless steel bow pull

Compressor:

Powder-coated steel

Manifold:

Powder-coated steel

Tubing:

Soft PVC

Compressor tube/air filter:

Soft PVC tube and polypropylene filter

Connectors (vary by type):

HDPE, nylon (white), polypropylene (white) and polycarbonate (clear)