


Ventaire™ Endoscopy System

Scope Drying Unit



Healthcare-associated infections (HAIs) affect one out of every 31 hospital patients at any given time, according to the Department of Health and Human Services' Agency for Healthcare Research and Quality. Over one million HAIs occur every year in the US, leading to loss of life and rising healthcare costs. A drying device that forces air through the interior lumens of the endoscope can reduce the level of moisture inside the scope and, therefore, the chance of contamination and infection.





Safe and thorough drying system

The Ventaire Endoscopy System includes products that assist endoscopy staff in the safe and thorough drying, tracking, and storage of reprocessed scopes, including scope drying and tracking cabinets and scope drying units. Both products provide consistent and continuous delivery of HEPA-filtered, pressurized air into endoscope channels to keep scopes dry.

Flexible, compact, and cost effective

The scope drying unit provides a flexible and cost-effective solution for thoroughly drying the interior lumens of endoscopes without the need for a dedicated drying cabinet. The compact size of the scope drying unit contributes to its flexibility – it can be placed on a counter or other surface, including carts, or mounted on a wall, IV pole, or scope pole. Scope drying units are available in two versions: one version has a compressor that supplies air and the other version connects to a building's air supply. The scope drying unit's compressor is balanced for low-vibration operation.

Once the unit's drying cycle is complete, the endoscope can be placed in a scope drying and tracking cabinet for continued drying and venting or stored in an endoscope cabinet. The scope drying unit works with all major endoscope brands.



The scope drying unit collects scope drying information, including date and time of operation, and identifies the operator and scope. Settings can be customized to comply with an organization's protocols and workflow. The display notifies operators when tubes and filters need replacing



The Ventaire Endoscopy System is a family of products that facilitate the storage, drying, and transport of endoscopes and includes drying and tracking cabinets, portable and attached scope drying units, and transport carts




Flexible in application, scope drying units can be placed on a counter or other surface, including Pace carts, or mounted on a wall, IV pole, or scope pole



The scope drying unit's barcode scanner and an easy-to-read, easy-to-use interface streamline operation

Scope drying units come in two versions: One has self-contained air supply, and the other connects to a building's air supply

Two scopes can be dried simultaneously



Connect to the InnerSpace Cloud platform and monitor scopes remotely, receive notifications, customize reports, and upload and download lists of users and scopes

SVDU2

Exterior: 7.5" w x 8" d x 10" h

Display screen: 5" diagonal length

Weight with compressor: 14.3 lb

Weight without compressor: 8.1 lb

Power cord length: 10 ft

Decibels (with compressor): 71 dB



Scope Drying Unit

The scope drying unit is available with or without connection to a building's air supply. Self-contained units have a compressor; units that connect to a building's air supply include a quick-disconnect plug.

SVDU2 – Scope Drying Unit, 2 drying ports, with compressor

SVDU2-H – Scope Drying Unit, 2 drying ports, without compressor (uses building's air supply).
Not to exceed 10 PSI

Accessories

A universal mounting bracket comes with a scope drying unit. The bracket can be used to mount the unit to a wall or an IV or scope pole attached to a Pace cart. Order IV or scope pole separately.

SIVPP – IV Pole for Pace cart

SIVPP-UAR – IV Pole for Pace cart, attaches to existing accessory rail post

SSCPP – Scope Pole for Pace cart, holds 2 scopes

SSCPP-UAR – Scope Pole for Pace cart, holds 2 scopes, attaches to existing accessory rail post

Replacement products

Tubes

Manifold and manifold auxiliary water tubes connect to the scope drying unit to effectively vent HEPA-filtered air through endoscope channels. Manifold tubes are for endoscopes without water jet ports; manifold auxiliary water tubes are for endoscopes with water jet ports. Tubes should be replaced every 24 hours.

SVDMT1 - Manifold Tube, package of 1

SVDWAT1 - Manifold Auxiliary Water Tube, package of 1

Connector Kits

Single-use connectors attach an endoscope's suction and air water valve ports to the scope drying unit to dry scopes. The easy-to-attach connectors support infection prevention protocols and are compatible with major endoscope brands. Kits include 100 connectors.

To determine which connector kit is compatible with your endoscope brands, go to

<https://innerspacehealthcare.com/Ventaire>.

HEPA Filter

Replace the filter every six months.

SVDUAF – Replacement HEPA Air Filter



InnerSpace specializes in solutions for all areas and departments within a healthcare facility.

We offer carts, cabinets, open storage, and accessories for supply storage, products with embedded technology for improved efficiency and security, and the Ventaire™ Endoscopy System for the storage, drying, transport, and reprocessing of endoscopes.

InnerSpace is here to help you specify the best products for your needs, whether a standard option or customized solution. We offer complimentary storage analysis and design services. Contact our Customer Service team with your requests or to connect with a sales representative, and find InnerSpace products on GSA Advantage, Premier, and Vizient contracts.



1239 Comstock Street | Marne, Michigan 49435
888-435-2256 phone | 616-743-5944 fax
innerspacehealthcare.com

© 2024 Solaire Medical, Marne, MI 49435
® InnerSpace is a registered trademark of Solaire Medical LLC, Marne, MI 49435
LIT-DRYINGUNIT-BR