
Accessory Product Catalog

Fine-tune carts and cabinets with the accessories you need, located where you need them.

Attached Accessories..... page 2

Accessories attached to the exteriors of carts keep essential tools within easy reach, saving time and providing easy access to essential supplies and equipment.

Procedure-specific Accessories page 15

Accessories target specific procedures for departments such as endoscopy and cardiac cath labs.

Interior Storage page 18

A modular system of shelves, dividers, trays, and baskets lets you arrange the interior configuration of your cabinets and carts for different types and sizes of supplies and easily reconfigure as supplies and needs change.

Locks..... page 27

Select the type of security you prefer with keyed and keyless options, including proximity locks.

Casework Finishing Accessories page 29

Add countertops and finishing touches to the installation of Evolve casework applications.

Side-Mounted Accessory Rail (SAMRSP, SAMRSUP, SAMRF24)

■ PACE ■ ROAM



Accessory rail attaches to the side of a Pace and Roam cart to hold a variety of accessories. Attachment hardware included.

Notes

Three rail styles work with different accessories and products. A notched rail accepts side-mounted accessories on a Pace cart. A universal rail accepts customers' own accessories on a Pace cart. A fixed rail accepts accessories on a Roam cart.

- Specify SAMRSP for notched adjustable side-mounted rail on Pace cart
- Specify SAMRSUP for universal adjustable side-mounted rail on Pace cart
- Specify SAMRF24 for fixed rail on Roam cart

Dimensions

- Pace cart: 17.75"w x 1.50"h
- Roam cart: 23.50"w x 1.50"h

Materials

- Notched rail: ABS thermoplastic polymer
- Universal rail: powder-coated steel
- Fixed rail: clear anodized aluminum

Weight capacity

- 25 lbs

Back-Mounted Accessory Rail (SAMRBP, SAMRBPN)

■ PACE



Accessory rail attaches to the back of a Pace cart to hold a variety of accessories. Attachment hardware included.

Notes

- Specify SAMRBP for adjustable back-mounted rail for Pace cart
- Specify SAMRBPN for adjustable back-mounted rail for narrow Pace cart

Dimensions

- Pace cart: 26.25"w x 1.50"h
- Narrow Pace cart: 14.50"w x 1.50"h

Materials

- Powder-coated steel

Weight capacity

- 20 lbs

Top-Mounted Accessory Rail (SPUAR1, SPUAR2, SPNUAR1, SPNUAR2)

PACE



Accessory rail attaches to the back rails on a Pace cart to hold a variety of accessories. Single rail has 1 rail with posts; double rail has 2 rails with posts. Maximum rail height is 24" above the cart surface. Attachment hardware included.

Notes

- Specify SPUAR1 for accessory rail with 1 rail
- Specify SPUAR2 for accessory rail with 2 rails
- Specify SPNUAR1 for accessory rail with 1 rail for narrow Pace cart
- Specify SPNUAR2 for accessory rail with 2 rails for narrow Pace cart
- Specify SPUARAIL for additional accessory rail
- Specify SPNUARAIL for additional accessory rail for narrow Pace cart
- Included with Pace Anesthesia, Difficult Airway, and IV Prep Carts

Dimensions

- Pace cart: 27.50"w x 4.7"h (per rail)
- Narrow Pace cart: 15.65"w x 4.7"h (per rail)

Materials

- Powder-coated steel

Weight capacity

- 20 lbs

Adjustable Work Surface (SPAWS, SPNAWS)

PACE



Pneumatic work surface attaches to the back rails on a Pace cart and provides a 12" height-adjustment range. Attachment hardware and back rails included.

Notes

- Specify SPAWS for Pace cart
- Specify SPNAWS for narrow Pace cart

Dimensions

- 17.75"w x 21"d (15.86"w x 19.61"d usable space)

Materials

- ABS thermoplastic polymer surface, aluminum extrusion

Weight capacity

- 30 lbs

Breakaway Lock Bar (SPLB, SPSBL)

■ PACE



Lock bar attaches to a Pace cart to provide additional security of cart contents. It can be sized to the height of a Pace cart or single drawer. Attachment hardware included.

Notes

- Specify SPLB for full-height lock
- Specify SPSBL for single-drawer lock
- Included with Pace Code Cart

Dimensions

- Full-height lock: 2"w (height will vary depending on height of cart)
- Single-drawer lock: 2"w x 3"h

Materials

- Powder-coated steel

Catheter Holder (SPCH)

■ PACE



Catheter holder attaches to the side rails on a Pace cart. Holds up to 20 pounds. Attachment hardware and side rails included.

Dimensions

- 8"w x 4"d x 19"h
- 4" space between fixed divider

Materials

- Powder-coated steel

Chart Holder (SCHP, SCH24)

■ PACE ■ ROAM



Chart holder attaches to the side rail on a Pace cart or to a fixed side rail on a Roam cart. Holds up to 5 pounds. Attachment hardware and side rail included.

Notes

- Specify SCHP for Pace cart
- Specify SCH24 for Roam cart
- Included with Pace Pharmacy Cart

Dimensions

- 13.43"w x 3.25"d x 7.75"h

Materials

- Powder-coated steel

Cord Wrap (SCWP, SCWPN)

■ PACE



Cord wrap attaches to the back rail on a Pace cart. Attachment hardware and back rail included.

Notes

- Specify SCWP for Pace carts
- Specify SCWPN for narrow Pace cart
- Order power strip (SPSP) separately

Dimensions

- 12"w x 2.5"d x .75"h

Materials

- HDPE cord wrap

Defibrillator Shelf (SDSP)

■ PACE



Defibrillator shelf attaches to the back of a Pace cart. It rotates 360 degrees (in 45-degree increments) and holds up to 25 pounds. Attachment hardware included.

Notes

- Cannot be attached to narrow Pace cart
- Included with Pace Code Cart

Dimensions

- 17"w x 13.75"d
- 14" clearance between top of cart and bottom of shelf

Materials

- HDPE shelf, powder-coated steel pole and brackets

Disinfectant Holder (SRDHP, SRDH24)

■ PACE ■ ROAM



Round disinfectant holder attaches to the side rail on a Pace cart or to a fixed side rail on a Roam cart. Holds containers up to 6" in diameter. Attachment hardware and side rail included.

Notes

- Specify SRDHP for Pace cart
- Specify SRDH24 for Roam cart
- Included with Pace Isolation and Pharmacy Carts

Dimensions

- 6.37"w x 4"d x 6"h

Materials

- Powder-coated steel with Velcro strap

Embolectomy Storage Rack (SESR)

■ ROAM



Storage rack attaches to the fixed side rails on a Roam cart to organize catheters. Holds up to 25 pounds. 5 fixed dividers, attachment hardware, and side rails included.

Dimensions

- 18"w x 3.5"d x 30"h
- 2.93" between dividers

Materials

- Powder-coated steel

Flip-Up Work Surface (SFWSP, SFWS24)

■ PACE ■ ROAM



Work surface attaches to the side rails on a Pace cart or to the fixed side rails on a Roam cart. Surface has seamless construction for ease of cleaning. Holds up to 30 pounds. Attachment hardware and side rails included.

Notes

- Cannot be used with narrow Pace cart
- Specify SFWSP for Pace cart
- Specify SFWS24 for Roam cart
- Included with Pace Pharmacy Cart

Dimensions

- 17"w x 13.75"d

Materials

- ABS thermoplastic polymer work surface

Glove Box Holder (S1GBHP, S1GBH24, S3GBHP, S3GBH24)

■ PACE ■ ROAM



Glove box holder attaches to the side rails on a Pace cart or to a fixed side rail on a Roam cart. Comes in 2 sizes and holds up to 5 pounds. Attachment hardware and side rails included.

Notes

- Specify S1GBHP for single box holder for Pace cart
- Specify S1GBH24 for single box holder for Roam cart
- Specify S3GBHP for triple glove box holder for Pace cart
- Specify S3GBH24 for triple glove box holder for Roam cart
- Included with Pace Isolation, IV Prep, and Pharmacy Carts

Dimensions

- Single box holder: 11.5"w x 3.6"d x 4.5"h
- Triple box holder: 11.5"w x 3.6"d x 14.5"h

Materials

- Powder-coated steel

IV Pole (SIVPP)

■ PACE



Adjustable IV pole attaches to the back of a Pace cart to hold and secure 2 IV bags. Pole has a 22" height-adjustment range and holds up to 10 pounds. Attachment hardware included.

Notes

- IV pole mounts anywhere along back channel of cart. Mounted at its highest position, the pole adjusts from 25" to 47" above cart top
- Included with Pace Code, Difficult Airway, and IV Prep Carts

Materials

- Stainless steel pole and steel tube base

LED Task Light (SPTL, SPNTL)

■ PACE



Adjustable LED task light attaches to the back rail on a Pace cart. It has a 26"-long flexible neck. Attachment hardware and back rail included.

Notes

- Specify SPTL for Pace cart
- Specify SPNTL for narrow Pace cart

Materials

- Powder-coated steel bracket

Monitor Arm (SFPMAP)

■ PACE



Monitor arm attaches to the upper rail on a Pace cart to hold a flat panel monitor up to 20 pounds and 27" wide. It has full tilt motion and a 25.75" height-adjustment range. Attachment hardware included.

Notes

- Cannot be used with narrow Pace cart
- Order accessory rail (SPUAR1) separately

Materials

- Powder-coated steel mount and aluminum, steel, and plastic arm

Multi-Roll Tape Dispenser (SSTD)

■ PACE



Freestanding tape dispenser has four individual sections and holds multiple rolls of tape with a total dispensing width of 9.5".

Notes

- Included with Pace Anesthesia Cart

Materials

- Powder-coated steel

Oxygen Tank Holder (SOTHP, SOTH24)

■ PACE ■ ROAM



Oxygen tank holder attaches to the side rails on a Pace cart or to the fixed side rails on a Roam cart and holds tanks up to 4" in diameter. Accommodates A, B, C, D, and E type tanks up to 12 pounds. Attachment hardware and side rails included.

Notes

- Specify SOTHP for Pace cart
- Specify SOTH24 for Roam cart
- Included with Pace Code Cart

Dimensions

- 6"w x 5"d x 19.50"h

Materials

- Powder-coated steel

Power Strip (SPSP, SPSEUP, SPSUKP)

■ PACE



Power strip attaches to the back rail on a Pace cart. The 120-volt medical-grade power strip has 6 outlets, on/off switch, and 15' cord. Attachment hardware and back rail included.

Notes

- Cannot be used with narrow Pace cart
- Specify SPSP for power strip
- Specify SPSEUP for power strip for European Union
- Specify SPSUKP for power strip for United Kingdom
- Order cord wrap (SCWP) separately

Dimensions

- 13.75"w x 1.75"d x 2.5"h

Pull-Out Keyboard Tray (SPKBT)

■ PACE



Keyboard tray fits into the top of a Pace cart in place of a 9"-high drawer to hold standard-sized keyboard. It has a 6" height-adjustment range and holds up to 10 pounds. Mouse pad and attachment hardware included.

Notes

- Cannot be used in narrow Pace cart

Dimensions

- 21"w x 10.5"d

Pull-Out Work Surface (SPPOWS, SPNPOWS)

■ PACE



Work surface fits into the top of a Pace cart to provide additional work space. Surface pulls out on the left or right side of a Pace cart and on the left side of a narrow Pace cart.

Notes

- Specify SPPOWS for Pace cart
- Specify SPNPOWS for narrow Pace cart

Dimensions

- Pace cart: 12"w x 10.30"d
- Narrow Pace cart: 9.8"w x 10.30"d

Materials

- HDPE

Weight capacity

- 10 lbs

Push Handle (SPAPH)

■ PACE



Push handle attaches to the side of a Pace cart. It is designed to push/pull from the side or top of the handle. Attachment hardware included.

Notes

- Included with Pace Anesthesia, Code, Difficult Airway, and Pharmacy Carts

Dimensions

- 16.25"w x 3.6"d x 9.25"h
- 3.5" interior depth

Materials

- Stainless steel

Replacement Top for Pace Cart (SPRT, SPNRT)

■ PACE



Removable and antimicrobial replacement top surface for a Pace cart.

Notes

- Specify SPRT for Pace cart
- Specify SPNRT for narrow Pace cart

Dimensions

- Pace cart: 27.21"w x 18.5"d
- Narrow Pace cart: 15.44"w x 18.5"d

Resuscitation Board (SRBP)

■ PACE



Resuscitation board attaches to the back of a Pace cart. Attachment hardware included.

Notes

- Cannot be used with narrow Pace cart
- Included with Pace Code Cart

Dimensions

- 20"w x 15"h

Materials

- HDPE board, powder-coated steel brackets

Sharps Container (SSCP, SSC24)

■ PACE ■ ROAM



Locking 5-quart Sharps container attaches to a side rail on a Pace cart or to a fixed side rail on a Roam cart. Holds up to 25 pounds. Attachment hardware and side rail included.

Notes

- Specify SSCP for Pace cart
- Specify SSC24 for Roam cart
- Included with Pace Anesthesia, Difficult Airway, and IV Prep Carts

Dimensions

- 10.5"w x 5"d x 11.75"h

Materials

- Plastic container, powder-coated steel holder

Side Storage Bin (SSSBP, SSSB24)

■ PACE ■ ROAM


Storage bin attaches to the side rail on a Pace cart or to a fixed side rail on a Roam cart. Holds up to 5 pounds. Attachment hardware and side rail included.

Notes

- Specify SSSBP for Pace cart
- Specify SSSB24 for Roam cart

Dimensions

- 9.5"w x 5.5"d x 6"h

Materials

- Powder-coated steel

Sloped Top for Roam Cart (SR1GST, SR2GST, SR3GST, SR4GST)

■ ROAM


Sloped top attaches to a Roam cart with hinged door. Attachment hardware included.

Notes

- Specify SR1GST for Roam 1 cart
- Specify SR2GST for Roam 2 cart
- Specify SR3GST for Roam 3 cart
- Specify SR4GST for Roam 4 cart
- Integrated sloped tops standard on Roam carts with roll-top doors

Dimensions

- Roam 1: 19.59"w x 28.23"d x 5"h
- Roam 2: 36.81"w x 28.23"d x 5"h
- Roam 3: 56.43"w x 28.23"d x 5"h
- Roam 4: 78.72"w x 28.23"d x 5"h

Materials

- Powder-coated steel

SneezeGuard Shield (SPSGS)

■ PACE


SneezeGuard shield attaches to the back of a Pace cart to provide a layer of protection from airborne transmission. Shield adjusts 360° and has a 21.5" height-adjustment range. Attachment hardware included.

Dimensions

- 20"w x 18.75"h x .25" thick

Materials

- Acrylic plastic shield, powder-coated steel pole

Suction Pump Holder (SSPHP, SSPH24)

■ PACE ■ ROAM



Suction pump holder attaches to the side of a Pace cart or to a fixed side rail on a Roam cart. Holds up to 20 pounds. Attachment hardware included.

Notes

- Specify SSPHP for Pace cart
- Specify SSPH24 for Roam cart
- Included with Pace Code Cart

Dimensions

- 19"w x 6.75"d x 10"h

Materials

- Powder-coated steel

Surface Rail (SPWAR, SPNWAR)

■ PACE



Surface rail attaches to the sides of a Pace cart to keep supplies and equipment from falling off surface. Attachment hardware included.

Notes

- Specify SPWAR for use with Pace cart
- Specify SPNWAR for use with Narrow Pace cart
- Included with Pace Anesthesia Cart

Dimensions

- Pace cart: 29"w x 18"d X 4"h
- Narrow Pace cart: 17"w x 18"d x 4"h
- 2" clearance between top of cart and bottom of rail

Materials

- Stainless steel

Tilt Bin Module (STB1, STB4, STB5, STB6)

■ PACE



Modules come in 1-, 4-, 5-, or 6-bin configurations and attach to the upper rail on a Pace cart. Bin module is enclosed with clear plastic fronts. Attachment hardware included.

Notes

- Specify STB1 for 1-compartment tilt bin, for use with narrow Pace cart only
- Specify STB4 for 4-compartment tilt bin
- Specify STB5 for 5-compartment tilt bin
- Specify STB6 for 6-compartment tilt bin
- 4-, 5-, and 6-bin modules cannot be used with narrow Pace cart
- Order accessory rail (SPUAR1, SPUAR2, SPNUAR1, SPNUAR2) separately
- 4- and 5-bin module included with Pace Anesthesia Cart; 4-bin module included with Pace Difficult Airway and IV Prep Carts

Dimensions

- STB1: 5.5"d x 5"w x 6.5"h
- STB4: 6.5"d x 23.5"w x 8"h
- STB5: 5.5"d x 23.5"w x 6.5"h
- STB6: 3.5"d x 23.5"w x 4.5"h

Materials

- ABS thermoplastic polymer body, acrylic thermoplastic polymer bin

Upper Utility Shelf (SPUUS, SPUUS-P, SPNUUS, SPNUUS-P)

■ PACE



Utility shelf attaches to upper posts on a Pace cart and holds up to 20 pounds. Attachment hardware included.

Notes

- Specify SPUUS for shelf only for Pace cart
- Specify SPUUS-P for shelf with posts for Pace cart
- Specify SPNUUS for shelf only for narrow Pace cart
- Specify SPNUUS-P for shelf with posts for narrow Pace cart
- Included with Pace Cast Cart

Dimensions

- Shelf for narrow Pace cart: 14.75"w x 12"d x 4.5"h
- Shelf for Pace cart: 26.75"w x 12"d x 4.5"h

Materials

- Powder-coated steel

Waste Basket (SWBP, SWB24)

PACE ROAM



Waste basket attaches to a side rail on a Pace cart or a fixed side rail on a Roam cart and holds up to 15 pounds. Attachment hardware and side rail included.

Notes

- Specify SWBP for Pace carts
- Specify SWB24 for Roam cart
- Included with Pace Anesthesia, Cast, and Difficult Airway Carts

Dimensions

- 11.5"w x 9"d x 13.25"h

Materials

- Plastic basket and powder-coated steel holder

Wire Basket (SSMWPB, SSMWB24)

PACE ROAM



Wire basket attaches to the side rail on a Pace cart or to a fixed side rail on a Roam cart and holds up to 5 pounds. Attachment hardware and side rail included.

Notes

- Specify SSMWPB for Pace cart
- Specify SSMWB24 for Roam cart

Dimensions

- 12.25"w x 6"d x 6"h

Materials

- Wire with powder-coated steel

Cath Manager® (S7CM, S9CM)

■ ROAM ■ EVOLVE



Attaches to a Metal Shelf for Cath Managers or a Plastic Shelf in a Roam cart or Evolve cabinet with FlexCell interior to manage catheters. Each inline hook holds 10-15 catheter packages and can be adjusted horizontally along shelf. Uses 5 FlexCells. Product labeling and attachment hardware included.

Notes

- Specify S7CM for 7-hook cath manager; attaches to 1- or 2-compartment 19.25"d shelf
- Specify S9CM for 9-hook cath manager; attaches to 1- or 2-compartment 27.25"d shelf
- Specify S7CMCARDS for additional 7-hook cath manager cards (package of 25)
- Specify S9CMCARDS for additional 9-hook cath manager cards (package of 25)

Cath Manager Module (S19M1CM4, S19M2CM7, S27P1CM3, S27M2CM7)

■ ROAM ■ EVOLVE



Module includes a shelf and cath managers for use in a Roam cart or Evolve cabinet with FlexCell interior. Uses 5 FlexCells.

Notes

- Specify S19M1CM4 for 4 7-hook cath manager module; attaches to 1-compartment 19.25"d metal shelf
- Specify S19M2CM7 for 7 7-hook cath manager module; attaches to 2-compartment 19.25"d metal shelf
- Specify S27P1CM3 for 3 9-hook cath manager module; attaches to 1-compartment 27.25"d plastic shelf
- Specify S27M2CM7 for 7 9-hook cath manager module; attaches to 2-compartment 27.25"d metal shelf
- Specify S7CMCARDS for additional 7-hook cath manager cards (package of 25)
- Specify S9CMCARDS for additional 9-hook cath manager cards (package of 25)

Fan Ventilation Unit (SFVU)

■ ROAM ■ EVOLVE



For use in Evolve scope cabinets and Roam scope carts. Includes HEPA filter and 10-amp, 125VAC 10-foot-long power cord. Attachment hardware included.

Notes

- Order replacement filters separately (SHEPFIL for 1 filter or SHEPFIL3 for 3 filters)

Specifications

- Voltage Rating: 115VAC
- Watts: 13.0W
- Amps: 0.19A

Dimensions

- 7"w x 5"d x 5"h

Hang-On Hangers (SH0)

Use hangers for hanging packages that are not pre-punched or are torn. Each hanger set includes 9 sheets with 20 hangers per sheet. Hangers have self-adhesive backing.

Pegboard Hooks



Hooks attach to a pegboard or a top-mounted accessory rail on Pace carts to hold and organize supplies.

Notes

- Specify SRPB4H for 4" hook
- Specify SRPB6H for 6" hook
- Specify SRPB8H for 8" hook
- Specify SRPB11H for 11" hook

Scope Manager® (SRSH1, SRSH2)



For use in Evolve scope cabinets and Roam scope carts. Metal bracket and attachment hardware included.

Notes

- Specify SRSH1 for single scope manager
- Specify SRSH2 for double scope manager

Materials

- ABS thermoplastic polymer with powder-coated metal bracket

Side-Mounted Scope Cabinet (SSMSCP)

■ PACE



Scope cabinet attaches to the side rails on a Pace cart to hold scopes up to 39" in length. Cabinet has a vented back, hinged polycarbonate door with keyed lock, and includes 3 scope managers with cord retainers. Attachment hardware and side rails included.

Notes

- Cannot be used with narrow Pace carts or Pace 21 cart
- Included with Pace Difficult Airway Cart

Dimensions

- 14.5"w x 12"d x 48"h

Materials

- HDPE top and bottom
- Extruded aluminum corner posts
- Alumatek panels
- Polycarbonate door
- Powder-coated steel attachment brackets

Slat Wall Hooks



Hooks attach to a slat wall to hold and organize supplies.

Notes

- Specify SRSW4H for 4" hook
- Specify SRSW6H for 6" hook
- Specify SRSW8H for 8" hook
- Specify SRSW10H for 10" hook
- Specify SRSW12H for 12" hook

Suture Module (SPSM)

■ ROAM ■ EVOLVE



Module includes 6 pull-out shelves to store and access suture boxes and other small supplies in a Roam cart or Evolve cabinet with FlexCell interior. Each shelf holds 8 2.25"-deep suture boxes for a total of 48 boxes. Uses 18 FlexCells.

Notes

- Order suture shelf dividers separately (SSSHDI)
- Order additional shelves separately (SSSH)
- Suture module holds up to 12 shelves (6 per side)

Dimensions

- 14"w x 23.5"d x 26.5"h

Materials

- Powder-coated steel on top and bottom and clear anodized aluminum extrusions and slatwall

6"H Basket, Clear (SB6)

PACE ROAM EVOLVE



For use in a Roam cart, Evolve cabinet, or Pace cart with FlexCell interior. Perforated bottoms and sides keep dust out of baskets. Can be easily removed to clean cart or cabinet interior. Integrated safety stops prevent basket from coming out of cart or cabinet. Holds up to 40 pounds. Uses 4 FlexCells.

Notes

- Baskets can be organized using adjustable dividers; order separately:
 - Specify SSD6 for clear short adjustable divider
 - Specify SLD6 for clear long adjustable divider

Dimensions

- Exterior: 15.75"w x 23.62"d x 5.20"h
- Interior: 13.86"w x 21.73"d x 5.03"h

9"H Basket, Clear (SB9)

PACE ROAM EVOLVE



For use in a Roam cart, Evolve cabinet, or Pace cart with FlexCell interior. Perforated bottoms and sides keep dust out of baskets. Can be easily removed to clean cart or cabinet interior. Integrated safety stops prevent basket from coming out of cart or cabinet. Holds up to 40 pounds. Uses 6 FlexCells.

Notes

- Baskets can be organized using adjustable dividers; order separately:
 - Specify SSD9 for clear short adjustable divider
 - Specify SLD9 for clear long adjustable divider

Dimensions

- Exterior: 15.75"w x 23.62"d x 8.32"h
- Interior: 13.43"w x 21.31"d x 8.15"h

Label Holder (SLH35, SLH25)




Label holder for 2" x 3.5" cards clips to the front of all trays and baskets. Label holder for 1.5" x 2" cards clips to an interior tray or basket divider. Both come in packages of 10. Cards not included.

Notes

- Specify SLH25 for divider label holder for 1.5" x 2" cards
- Specify SLH35 for label holder for 2" x 3.5" cards

Snap-Off Divider Kit (SSODK)



For use in all 3"-high trays. Clear snap-off divider is adjustable and held in place by a divider clip that attaches anywhere along a tray divider and can hold 2 snap-off dividers (1 per side). Kit includes 5 snap-off dividers and 5 divider clips.

Dimensions

- Snap-off divider: 5.21"w x 2.43"h, snaps off in .63" increments
 - Total of 7 snap-off sections
- Divider clip: 1.13"w x 1"d x 2"h

3"H Tray (SCT3, ST3, STD3)

■ PACE
 ■ ROAM
 ■ EVOLVE



For use in a Roam cart, Evolve cabinet, or Pace cart with FlexCell interior or in a 3"-high Pace drawer. Tray has seamless construction for ease of cleaning. Holds up to 40 pounds. Uses 2 FlexCells.

Notes

- Specify SCT3 for clear pharmacy tray
- Specify ST3 for grey tray
- Pharmacy tray included with Pace Pharmacy Cart

Trays can be organized using adjustable tray dividers; order separately:


- Specify STD3 for grey preconfigured tray with 3 short and 3 long dividers
- Specify SSD3 clear short tray divider
- Specify SLD3 clear long tray divider
- Specify SSOD for adjustable snap-off divider and SSODC for a snap-off divider clip

Dimensions

- Exterior: 15.75"w x 23.62"d x 2.63"h
- Interior: 13.86"w x 21.73"d x 2.47"h

6”H Tray, Grey (ST6, STD6, SST6)

PACE ROAM EVOLVE



For use in a Roam cart, Evolve cabinet, or Pace cart with FlexCell interior or in a 6"-high Pace drawer. Tray has seamless construction for ease of cleaning. Holds up to 40 pounds. Uses 4 FlexCells.

Notes

- Specify ST6 for grey tray
- Specify SST6 for smooth tray and SSTL6 for smooth tray lid

Trays can be organized using adjustable tray dividers; order separately:

- Specify STD6 for preconfigured tray with 3 short and 3 long dividers
- Specify SSD6 for clear short tray divider
- Specify SLD6 for clear long tray divider

Dimensions

- Exterior: 15.75"w x 23.62"d x 5.20"h
- Interior: 13.86"w x 21.73"d x 5.03"h

9"H Tray, Grey (ST9, STD9)

PACE ROAM EVOLVE



For use in a Roam cart, Evolve cabinet, or Pace cart with FlexCell interior or in a 9"-high Pace drawer. Tray has seamless construction for ease of cleaning. Holds up to 40 pounds. Uses 6 FlexCells.

Notes

- Specify ST9 for grey tray

Trays can be organized using adjustable tray dividers; order separately:

- Specify STD9 for preconfigured tray with 1 short and 1 long divider
- Specify SSD9 for clear short tray divider
- Specify SLD9 for clear long tray divider

Dimensions

- Exterior: 15.75"w x 23.62"d x 8.32"h
- Interior: 13.43"w x 21.31"d x 8.15"h

3"H Narrow Tray, Grey (SNT3, SNTD3)

PACE EVOLVE



For use in a narrow 3"-high Pace drawer or Evolve FlexCell cabinet with offset column. Tray has seamless construction for ease of cleaning. Holds up to 40 pounds. Uses 2 FlexCells.

Notes

- Specify SNT3 for grey tray
- Trays can be organized using adjustable tray dividers; order separately:
 - Specify SNTD3 for preconfigured tray with 3 short and 1 long divider
 - Specify SSND3 for clear short tray divider
 - Specify SSD3 for clear long tray divider
 - Specify SSOD for adjustable snap-off divider and SSODC for snap-off divider clip

Dimensions

- Exterior: 11.81"w x 15.75"d x 2.63"h
- Interior: 9.80"w x 13.73"d x 2.47"h

6"H Narrow Tray, Grey (SNT6, SNTD6)

PACE EVOLVE



For use in a narrow 6"-high Pace drawer or Evolve FlexCell cabinet with offset column. Tray has seamless construction for ease of cleaning. Holds up to 40 pounds. Uses 4 FlexCells.

Notes

- Specify SNT6 for grey tray

Trays can be organized using adjustable tray dividers; order separately:

- Specify SNTD6 for preconfigured tray with 3 short and 1 long divider
- Specify SSND6 for clear short tray divider
- Specify SSD6 for clear long tray divider

Dimensions

- Exterior: 11.81"w x 15.75"d x 5.20"h
- Interior: 9.66"w x 13.60"d x 5.03"h

9"H Narrow Tray, Grey (SNT9, SNTD9)

PACE EVOLVE



For use in a narrow 9"-high Pace drawer or Evolve FlexCell cabinet with offset column. Tray has seamless construction for ease of cleaning. Holds up to 40 pounds. Uses 6 FlexCells.

Notes

- Specify SNT9 for grey tray

Trays can be organized using adjustable tray dividers; order separately:

- Specify SNTD9 for preconfigured tray with 1 short and 1 long divider
- Specify SSND9 for clear short tray divider
- Specify SSD9 for clear long tray divider

Dimensions

- Exterior: 11.81"w x 15.75"d x 8.32"h
- Interior: 9.50"w x 13.43"d x 8.15"h

Basic Metal Shelf (S19S1)

■ EVOLVE



Fits anywhere within a single compartment in a 19.25"-deep Evolve cabinet with FlexCell interior. Holds up to 25 pounds. Uses 1 FlexCell.

Dimensions

- 23.5"w x 16"d

Narrow Basic Metal Shelf (S19NS1)

■ EVOLVE



Fits in the narrow compartment of an Evolve cabinet with offset column. Holds up to 25 pounds. Uses 1 FlexCell.

Dimensions

- 11.5"w x 16"d

Metal Shelf for Cath Managers (S19M1, S19M2, S27M2)

■ ROAM ■ EVOLVE



Fits anywhere in a Roam cart or Evolve cabinet to hold cath managers. Holds up to 50 pounds. Uses 1 FlexCell.

Notes

- Specify S19M1 for 1-compartment 19.25"d cabinet
- Specify S19M2 for 2-compartment 19.25"d cabinet
- Specify S27M2 for 2-compartment 27.25"d cart or cabinet

Dimensions

- S19M1: 23.5"w x 16"d
- S19M2: 32.75"w x 16"d
- S27M2: 32.75"w x 23.75"d

Plastic Shelf (S27P1)

■ PACE ■ ROAM ■ EVOLVE

Fits anywhere in a Pace cart with FlexCell interior, Roam cart, or 27.25"-deep Evolve cabinet with 1 compartment. Can be used as a standalone shelf or to hold cath managers. Holds up to 50 pounds without cath manager; 25 pounds with cath manager. Uses 1 FlexCell.

Dimensions

- 15.75"w x 23.75"d

Adjustable Shelf Divider (S19DI, S27DI)

■ ROAM ■ EVOLVE



Adjustable and removable 5"-high plastic divider has 2 notches designed to hold a backstop for smaller boxes; order separately. Uses 4 FlexCells.

Notes

- Shelf backstops located 9.75" and 12.25" from front of shelf divider
- Specify S19DI for use with 19.25"d shelf
- Specify S27DI for use with 27.25"d shelf
- Specify S19BS1 for divider backstop for 1-compartment 19.25"d divider
- Specify S19BS2 for divider backstop for 2-compartment 19.25"d divider
- Specify S27BS1 for divider backstop for 1-compartment 27.25"d divider
- Specify S27BS2 for divider backstop for 2-compartment 27.25"d divider

Drawers (S19DR1, S19DR2, S27DR1, S27DR2)

■ ROAM ■ EVOLVE



For use with a Roam cart or Evolve cabinet with FlexCell interior to store and organize supplies. Includes keyed lock. Holds up to 40 pounds. Uses 5 FlexCells.

Notes

- Specify S19DR1 for 22.5"w x 6"h drawer for 1-compartment 19.25"d cabinet
- Specify S19DR2 for 31.75"w x 6"h drawer for 2-compartment 19.25"d cabinet
- Specify S27DR1 for 14.5"w x 6"h drawer for 1-compartment 27.25"d cart or cabinet
- Specify S27DR2 for 31.75"w x 6"h drawer for 2-compartment 27.25"d cart or cabinet

Dimensions

- S19DR1: interior, 18.76"w x 15.65"d x 6.19"h
- S19DR2: interior, 30.52"w x 15.65"d x 6.19"h
- S27DR1: interior, 13.32"w x 23.65"d x 6.19"h
- S27DR2: interior, 30.55"w x 23.65"d x 6.19"h

Materials

- Powder-coated steel

Boxed Catheter Insert (S19BC1, S19BC2, S27BC1, S27BC2, S27BC2-AT)

■ ROAM ■ EVOLVE



Boxed Catheter Insert with anti-tip mechanism shown

For use in a Roam cart or Evolve cabinet to hold catheters. Holds up to 50 pounds. Attachment hardware included. Uses 29 FlexCells.

Notes

- Specify S19BC1 for use in 1-compartment 19.25"d cabinet. Includes 3 adjustable dividers
- Specify S19BC2 for use in 2-compartment 19.25"d cabinet. Includes 5 adjustable dividers
- Specify S27BC1 for use in 1-compartment 27.25"d cart or cabinet. Includes 2 adjustable dividers
- Specify S27BC2 for use in 2-compartment 27.25"d cabinet. Includes 5 adjustable dividers
- Specify S27BC2-AT for use in 2-compartment 27.25"d cart. Includes anti-tip mechanism and 4 adjustable dividers (2 per section)

Dimensions

- S19BC1: exterior, 22"w x 16"d x 30-42.25"h;
interior, 20.48"w x 14.48"d x 26.24-38.24"h.
Extends 10.75" from front edge of cabinet
- S19BC2: exterior, 31"w x 16"d x 30-42.25"h;
interior, 29.48"w x 14.48"d x 26.24-38.24"h.
Extends 10.75" from front edge of cabinet
- S27BC1: exterior, 14.25"w x 23.50"d x 30-42.25"h;
interior, 12.73"w x 21.98"d x 26.24-38.24"h.
Extends 16.25" from front edge of cart or cabinet
- S27BC2: exterior, 31"w x 23.50"d x 30-42.25"h;
interior, 29.48"w x 21.98"d x 26.24-38.24"h.
Extends 16.25" from front edge of cabinet
- S27BC2-AT: exterior, 31"w x 21.50"d x 30-42.25"h;
interior, each compartment is 15.15"w x 21.40"d x 26.24-38.24"h.
Extends 12.75" from front edge of cart

Materials

- Insert for Roam cart: Powder-coated steel box and base, clear acrylic dividers on anti-tip mechanism
- Insert for Evolve cabinet: Melamine box, clear acrylic dividers

LED Lights (LED)

■ ROAM ■ EVOLVE

Lights attach to channels along the front of Roam carts and Evolve cabinets with hinged doors to provide interior lighting. Lights are cool white (6000K) and have an average life of 50,000 hours. Not available on Evolve upper cabinets with sloped tops.

Security Lock Box (SPSLB3, SPNSLB3)

■ PACE



For use with a 3"-high Pace drawer or 3"-high narrow Pace drawer. Security lock box is secured to the drawer and includes a keyed lock and recessed handle. Attachment hardware included.

Notes

- Specify SPSLB3 for Pace 3"h drawer
- Specify SPNSLB3 for Narrow Pace 3"h drawer
- 2 narrow security boxes can fit side-by-side in a 3"h Pace drawer
- Included with Pace Pharmacy Cart

Dimensions

- For 3"h drawer: 23.13"w x 15.53"d x 2.74"h
(interior: 22.73"w x 14.31"d x 1.31-2.74"h)
- For narrow 3"h drawer: 11.32"w x 15.53"d x 2.74"h
(interior: 10.92"w x 14.31"d x 1.31-2.74"h)


Materials

- Powder-coated steel



Keyed Lock (SLK)

PACE ROAM EVOLVE



For keyed access to a Roam cart, Pace cart, or Evolve cabinet. Includes 2 keys.

Notes

- Keyed lock standard with Pace carts

Keypad Lock (SLE)

PACE ROAM EVOLVE



For keypad access to a Roam cart, Pace cart, or Evolve cabinet.

Proximity Lock (SLP)

PACE ROAM EVOLVE



For badge access to a Roam cart, Pace cart, or Evolve cabinet for up to 2,000 users and administrators. Includes an LED low-battery indicator.

Notes

- Specify SSMART for pack of 5 smart cards for proximity locks

Keypad with Proximity Lock (SLD)

■ PACE ■ ROAM ■ EVOLVE



For badge and keypad access to a Roam cart, Pace cart, or Evolve cabinet for up to 2,000 users and administrators. Includes an LED low-battery indicator.

Notes

- Specify SSMART for pack of 5 smart cards for proximity locks

Countertop, Compression Board with Laminate (SEC19CT, SEC27CT)

■ EVOLVE



Laminated 1.125"-thick surface fits on an Evolve base cabinet. It has a 4"-high backsplash, finished ends, and Wilsonart Frosty White (1573) surface.

Notes

- Specify SEC19CT for 19.25"d cabinet; countertop is 19.75"d
- Specify SEC27CT for 27.25"d cabinet; countertop is 27.75"d
- Countertop can be specified to any length up to 96"
- Contact InnerSpace Customer Service to specify alternative laminate
- Optional 1.5"-thick countertop available in stainless steel and can be specified in any length up to 138"

Countertop Support Bracket (SEC19WSB, SEC27WSB, SES19WSB, SES27WSB)

■ EVOLVE

Bracket provides additional support for an Evolve countertop that extends beyond the base cabinet.

Notes

- Specify SEC19WSB for metal bracket used with laminate 19.25"d countertop
- Specify SEC27WSB for metal bracket used with laminate 27.25"d countertop
- Specify SES19WSB for stainless steel bracket used with stainless steel 19.25"d countertop
- Specify SES27WSB for stainless steel bracket used with stainless steel 27.25"d countertop

Finished End (TFE, UFE, BFE)

■ EVOLVE



Finishes the sides of an Evolve cabinet to match the cabinet's laminate.

Notes

- Specify TFE for infill panel used with tall cabinet
- Specify UFE for infill panel used with upper cabinet
- Specify BFE for infill panel used with base cabinet

Infill Panel (TIF, UIF, BIF)

■ EVOLVE

Fills the gaps between adjacent Evolve cabinets in a run. Panel height matches the height of the adjacent cabinet, and panel finish matches cabinets' laminate surfaces or brushed aluminum finishes. Maximum length is 96".

Notes

- Specify TIF for infill panel used with tall cabinet
- Specify UIF for infill panel used with upper cabinet
- Specify BIF for infill panel used with base cabinet

Sloped Top (THST, TRST, UHST, SHST, SRST)

■ EVOLVE



Sloped top attaches to an Evolve cabinet with a hinged or roll-top door. Attachment hardware included.

Notes

- 9.75"h sloped top for 27.25"d cabinet:
 - Specify 9THST for cabinet with hinged doors
 - Specify 9TRST for cabinet with roll-top door
- 7"h sloped top for 19.25"d cabinet:
 - Specify 7THST for cabinet with hinged doors
 - Specify 7TRST for cabinet with roll-top door
- 4.75"h sloped top for 13"d cabinet:
 - Specify 4UHST for cabinet with hinged doors
- 3"h vented sloped top for scope cabinet:
 - Specify 3SHST for cabinet with hinged doors
 - Specify 3SRST for cabinet with roll-top door

Materials

- Sloped top material and finish will match cabinet material and finish

