
# Evolve Scope Cabinet with Hinged Door

Evolve modular cabinets create a flexible storage system that organizes a hospital's most supply-intensive departments. A clean design and tight construction offer a refined aesthetic and support infection prevention efforts. The 19.25"-deep scope cabinet comes in 3 widths and 2 heights. It has a solidfront or tempered glass door, padded rear wall, removable drip tray, and adjustable scope managers and cord retainers. Material choices include compression board with a melamine or specified laminate surface, stainless steel, or nonporous AireCore with Glacier White finish, Brushed Aluminum finish, or specified laminate surface. Cabinet has a 4"-high toe kick.





## Evolve Scope Cabinet with Hinged Door



#### **Dimensions** (nominal)

Exterior (interior); scope hanging length: 26.75"w x 19.25"d x 60"h (25.25"w x 17"d x 54.5"h); 51" maximum 26.75"w x 19.25"d x 92"h (25.25"w x 17"d x 86.5"h); 83" maximum 36"w x 19.25"d x 92"h (34.5"w x 17"d x 86.5"h) (shown); 83" maximum 40"w x 19.25"d x 92"h (38.5"w x 17"d x 86.50"h); 83" maximum 60"h cabinet with 6"h toe kick is 62"h 92"h cabinet with 6"h toe kick is 94"h

#### **Interior Configuration**

60"h cabinet holds 5 scopes

92"h x 26.75"w cabinet holds 5 or 10 scopes

92"h x 36"w cabinet holds 8 or 16 scopes

92"h x 40"w cabinet holds 9 or 18 scopes

#### **Options**

- 92"h cabinet includes solid-front and glass door options; 60" cabinets have glass doors only
- Door can be hinged on right or left side of 26.75"w cabinet
- Fan Ventilation Unit
- · Keyed, InterConnect lock, and proximity lock upgrade for InterConnect lock
- Sloped top
- 6" toe kick
- InnerSpace Cloud for remote management and monitoring of cabinet security





### Product numbering system for Evolve Scope Cabinet with Hinged Door

Product Shown: SEC1936SGS

SEC	1 9	9 3 6	6 S	G		S	
Product Line Material	Depth	Width	Cabinet Type	Door Type	Hinge Location	Cabinet Configuration	
	26.75" Compre AireCor AireCor	ct Numbers w x 19.25"d x 60"h ession board: SEC1927 e: SEA1927GLB, SEA1 e and Brushed Aluminu ss steel: SES1927GLB,	927GRB um: SEB1927GLB,	. ,			
1	26.75" Compre	Product Numbers 26.75"w x 19.25"d x 92"h Compression board: SEC1927SGLS, SEC1927SGRS, SEC1927SSLS, SEC1927SSRS; with double scope holder: SEC1927SGLD, SEC1927SGRD, SEC1927SSLD, SEC1927SSRD					
	AireCor holder: AireCor scope h	AireCore: SEA1927SGLS, SEA1927SGRS (shown), SEA1927SSLS, SEA1927SSRS (shown); with double scope holder: SEA1927SGLD, SEA1927SGRD, SEA1927SSLD, SEA1927SSRD AireCore and Brushed Aluminum: SEB1927SGLS, SEB1927SGRS, SEB1927SSLS, SEB1927SSRS; with double scope holder: SEB1927SGLD, SEB1927SGRD, SEB1927SSLD, SEB1927SSRD					
	with do	Stainless steel: SES1927SGLS, SES1927SGRS, SES1927SSLS, SES1927SSRS; with double scope holder: SES1927SGLD, SES1927SGRD, SES1927SSLD, SES1927SSRD					
	36"w x Compre AireCor AireCor SEB193	Product Numbers36"w x 19.25"d x 92"hCompression board: SEC1936SGS, SEC1936SSS; with double scope holder: SEC1936SGD, SEC1936SSDAireCore: SEA1936SGS, SEA1936SSS (shown); with double scope holder: SEA1936SGD, SEA1936SSDAireCore and Brushed Aluminum: SEB1936SGS, SEB1936SSS; with double scope holder: SEB1936SGD, SEB1936SSDStainless steel: SES1936SGS, SES1936SSS; with double scope holder: SES1936SGD, SES1936SSD					
	Product 40"w x Compre AireCor AireCor	Product Numbers 40"w x 19.25"d x 92"h Compression board: SEC1940SGS, SEC1940SSS; with double scope holder: SEC1940SGD, SEC1940SSD AireCore: SEA1940SGS, SEA1940SSS (shown); with double scope holder: SEA1940SGD, SEA1940SSD AireCore and Brushed Aluminum: SEB1940SGS, SEB1940SSS; with double scope holder: SEB1940SGD,					
	SEB194 Stainles	40SSD ss steel: SES1940SGS,	SES1940SSS; wit	h double scope ho		D, SES1940SSD   innerspacehealthcare.com	

## Evolve Scope Cabinet with Hinged Door

#### Materials and finishes

#### Cabinet shell:

- Compression board: Thermally fused melamine over 0.75"-thick compression board or specified Wilsonart, Formica, or Pionite standard laminate
- AireCore: Aluminum composite with Glacier White finish or specified Wilsonart, Formica, or Pionite standard laminate
- · AireCore and Brushed Aluminum finish: Aluminum composite with Brushed Aluminum finish
- Stainless steel: 18-gauge 304 stainless steel with 18% chrome and 8% nickel analysis; exposed surfaces

#### Edge banding:

- Compression board and AireCore: White ABS
- · AireCore with Brushed Aluminum finish: Brushed Aluminum acrylic

#### Solid-front door:

- Compression board: Thermally fused melamine over 0.75"-thick compression board or specified Wilsonart, Formica, or Pionite standard laminate
- AireCore: Aluminum composite with Glacier White finish or specified Wilsonart, Formica, or Pionite standard laminate
- · AireCore and Brushed Aluminum finish: Aluminum composite with Brushed Aluminum finish
- Stainless steel: 18-gauge 304 stainless steel with 18% chrome and 8% nickel analysis

#### Glass-front door:

- Compression board: Thermally fused melamine over 0.75"-thick compression board or specified Wilsonart, Formica, or Pionite standard laminate with tempered glass and clear anodized aluminum trim
- AireCore: Aluminum composite with Glacier White finish or specified Wilsonart, Formica, or Pionite standard laminate with tempered glass and clear anodized aluminum trim
- AireCore and Brushed Aluminum finish: Aluminum composite with Brushed Aluminum finish with tempered glass and clear anodized aluminum trim
- Stainless steel: 18-gauge 304 stainless steel with 18% chrome and 8% nickel analysis and tempered glass

#### Hinges:

- Compression board, AireCore, and AireCore with Brushed Aluminum cabinet: Metal concealed European-style
- Stainless steel cabinet: Stainless steel

#### Door pulls:

- · Compression board, AireCore, and AireCore with Brushed Aluminum cabinet: Dull chrome-plated steel bow pull
- · Stainless steel cabinet: 5" stainless steel bow pull

