


Replacement Products

Connector Kit Example



Support scope drying and storage processes, safety protocols, and reduce the chances (and time and expense) of reprocessing expired endoscopes with the scope drying and tracking cabinet and scope drying unit. The drying system meets the 2021 Multisociety Task Force Scope Storage and Drying Guidelines.

Connector Kits

Single-use connectors attach an endoscope's suction and air water valve ports to the Ventaire system to effectively dry and vent scopes. The easy-to-attach connectors support infection prevention protocols and are compatible with major endoscope brands. The connectors have undergone ethylene oxide sterilization. Kits include 100 connectors.

To determine which connector kit is compatible with your endoscope brands, visit <https://innerspacehealthcare.com/Ventaire>.

Materials and Finishes

HDPE, white nylon, white polypropylene, clear polycarbonate