

# Ventaire™ Endoscopy System

## Replacement Products


Compressor Tube / Air Filter



HEPA Air Filter



HEPA Fan Filter



## Ventaire System HEPA Filters

The filter pushes HEPA filtered air through an endoscope. Compressor tube/air filters are used with scope drying and tracking cabinets; HEPA air filters are used with scope drying units. Replace filter every six months.

Scope drying and tracking cabinet  
Replacement compressor tube/air filter: SVCTAF  
(Package of 1)

Scope drying unit  
Replacement HEPA air filter SVDUAF  
(Package of 1)

## Materials and Finishes

Scope drying and tracking cabinet:

Compressor tube and air filter: Soft PVC tube and polypropylene filter

Scope drying unit:

Filter: Polypropylene

## HEPA Fan Filter

The HEPA filter is used in an InnerSpace® scope cabinet's fan unit. Replace the filter every six months.

Package of 1 (SHEPFIL)

Package of 3 (SHEPFIL3)