

# Scope Tracking Retrofit Unit

---



The scope tracking retrofit unit is installed in a previously purchased Evolve scope cabinet. It includes a scope tracking unit, attachment hardware, and replacement door with 10" touchscreen. The unit is designed for easy installation and start-up.

The scope tracking retrofit unit efficiently and cost-effectively manages endoscope storage. The clear and visual interface improves scope monitoring and reduces the chances (and time and expense) of reprocessing expired endoscopes. The tracking system identifies scope locations, storage times, and lengths of time before scopes expire. It also generates a report that records all scope activity within the cabinet.



- Support compliance by monitoring hanging times of clean scopes
- User friendly 10" touchscreen provides clear operation for monitoring endoscope status
- Monitor scope storage from touchscreen on cabinet or remotely via the InnerSpace Cloud Platform
- Modify or add to scope inventory at any time
- Customizable report captures scope information, including date and time, user, scope expiration times, and check-in and check-out log

# Scope Tracking Retrofit Unit

---

---

## Options

- Subscription to InnerSpace Cloud Services for remote monitoring of scopes

---

## Specifications


Electrical requirements for tracking unit:

- 120 VAC, 60Hz, 15-amp standard duplex receptacle
- Includes a 125V, 10-amp, 10-foot-long power cord

# Ventaire™ Endoscopy System

## Scope Tracking Retrofit Unit

### Product Number Reference



### Product Numbers

To order a tracking retrofit unit, first specify the replacement door size to match your cabinet size, and then specify the retrofit tracking unit.

#### Door for 26.75"w x 48"h left-hinged cabinet

- AireCore: SEA1927SGLWVD (glass door), SEA1927SSLWVD (solid door)
- AireCore with Brushed Aluminum: SEB1927SGLWVD (glass door), SEB1927SSLWVD (solid door)
- Stainless steel: SES1927SGLWVD (glass door), SES1927SSLWVD (solid door)

#### Door for 26.75"w x 48"h right-hinged cabinet

- AireCore: SEA1927SGRWVD (glass door), SEA1927SSRWVD (solid door)
- AireCore with Brushed Aluminum: SEB1927SGRWVD (glass door), SEB1927SSRWVD (solid door)
- Stainless steel: SES1927SGRWVD (glass door), SES1927SSRWVD (solid door)

#### Door for 26.75"w x 60"h left-hinged cabinet

- AireCore: SEA1927SGLBVD (glass door), SEA1927SSLBVD (solid door)
- AireCore with Brushed Aluminum: SEB1927SGLBVD (glass door), SEB1927SSLBVD (solid door)
- Stainless steel: SES1927SGLBVD (glass door), SES1927SSLBVD (solid door)

#### Door for 26.75"w x 60"h right-hinged cabinet

- AireCore: SEA1927SGRBVD (glass door), SEA1927SSRBVD (solid door)
- AireCore with Brushed Aluminum: SEB1927SGRBVD (glass door), SEB1927SSRBVD (solid door)
- Stainless steel: SES1927SGRBVD (glass door), SES1927SSRBVD (solid door)

#### Door for 26.75"w x 92"h left-hinged cabinet

- AireCore: SEA1927SGLVD (glass door), SEA1927SSLVD (solid door)
- AireCore with Brushed Aluminum: SEB1927SGLVD (glass door), SEB1927SSLVD (solid door)
- Stainless steel: SES1927SGLVD (glass door), SES1927SSLVD (solid door)

#### Door for 26.75"w x 92"h right-hinged cabinet

- AireCore: SEA1927SGRVD (glass door); SEA1927SSRVD (solid door)
- AireCore with Brushed Aluminum: SEB1927SGRVD (glass door); SEB1927SSRVD (solid door)
- Stainless steel: SES1927SGRVD (glass door); SES1927SSRVD (solid door)



# Ventaire™ Endoscopy System

## Scope Tracking Retrofit Unit

### Product Number Reference

S E		A	1	9	3	6	S	G	V D	
Product Line		Material	Depth		Width		Cabinet Type	Door Type	Hinge Location	Cabinet Configuration

### Product Numbers

To order a tracking retrofit unit, first specify the replacement door size to match your cabinet size, and then specify the retrofit tracking unit.

#### Left-hinged door for 36"w cabinet

- AireCore: SEA1936SGLVD (glass door); SEA1936SSLVD (solid door)
- AireCore with Brushed Aluminum: SEB1936SGLVD (glass door); SEB1936SSLVD (solid door)
- Stainless steel: SES1936SGLVD (glass door); SES1936SSLVD (solid door)

#### Left-hinged door for 40"w cabinet

- AireCore: SEA1940SGLVD (glass door); SEA1940SSLVD (solid door)
- AireCore with Brushed Aluminum: SEB1940SGLVD (glass door); SEB1940SSLVD (solid door)
- Stainless steel: SES1940SGLVD (glass door); SES1940SSLVD (solid door)

# Ventaire™ Endoscopy System

## Scope Tracking Retrofit Unit

### Product Number Reference



### Product Numbers

Scope Tracking Retrofit Unit – For 26.75" w left-hinged cabinet

- 5-scope unit: SVTURL5
- 10-scope unit: SVTURL10

Scope Tracking Retrofit Unit – For 26.75" w right-hinged cabinet

- 5-scope unit: SVTURR5
- 10-scope unit: SVTURR10

Scope Tracking Retrofit Unit – For 36" w cabinet

- 8-scope unit: SVTUR8
- 16-scope unit: SVTUR16

Scope Tracking Retrofit Unit – For 40" w cabinet

- 9-scope unit: SVTUR9
- 18-scope unit: SVTUR18