

Ventaire™ Endoscopy System

Replacement Products

Compressor Tube Replacement



Manifold Tube



Manifold Auxiliary Water Tube



Tubes

Compressor, manifold, and manifold auxiliary water tubes connect to components in the drying system of a scope drying and tracking cabinet or to the ports of a scope drying unit to effectively vent HEPA filtered air through endoscope channels. Manifold tubes are used with endoscopes that don't have water jet ports; manifold auxiliary water tubes are used with endoscopes with water jet ports. Tubes should be replaced every six months.

Tubes for scope drying and tracking cabinet

Compressor tube replacement (SVCTAF)

Includes air filter, package of 1

Triple manifold tube (SVTMT1)

Package of 1

Manifold tube (SVMT1)

Package of 1

Manifold auxiliary water tube (SVWAT1)

Package of 1

Tubes for scope drying unit

Triple manifold tube (SVTMT1)

Package of 1

Manifold tube (SVDMT1)

Package of 1

Manifold auxiliary water tube (SVDWAT1)

Package of 1

Materials and Finishes

Compressor tube and air filter:

Soft PVC tube and polypropylene filter

Manifold and manifold auxiliary water tubes:

Soft PVC