


ECP Cart with Drawers and 2 Center Columns (SRSDCECP)

Roam carts are lightweight and easy to maneuver. The FlexCell design and system of interior components adapt to the specific supply and procedure storage needs of departments throughout a healthcare facility. The Roam 3 ECP (extending cardiovascular procedure) Cart with Drawers and 2 Center Columns has a

stainless steel top that extends both left and right to increase surface space, 2 full-width drawers, and 8 FlexCells in each of the cart's 3 compartments. The drawers have full-length extension glides and integral pulls. The cart comes in a Soft White finish.



ECP Cart with Drawers and 2 Center Columns (SRSDCECP)



Materials and Finishes

Base structure: Powder-coated steel

Top structure: 304 Stainless Steel

Corner extrusions: Clear anodized aluminum

Corner bumpers: Thermoplastic polymer (ABS plastic)

FlexCell panels: Thermoplastic polymer (ABS plastic)

Drawers: Powder-coated steel

Casters: Nylon and steel

Dimensions (nominal) and Weight

Exterior: 57"w x 28.75"d x 36"h, 88.25"w when left and right pull-out surfaces are extended

Interior: 14.5" w (each compartment) x 23.75" d x 12.75"h

Weight: 330 lbs

Interior Configuration

24 FlexCells, 8 per compartment

2 full-width drawers

Options

Interior and side-mounted exterior accessories