


ECP Cart with Drawers and 2 Center Columns (SRSDCECP)

Roam carts are lightweight and easy to maneuver. The FlexCell design and system of interior components adapt to the specific supply and procedure storage needs of departments throughout a healthcare facility.

The Roam 3 ECP (extending cardiovascular procedure) Cart with Drawers and 2 Center Columns has a stainless steel top that extends both left and right to increase surface space, 2 full-width drawers, and 5 FlexCells in each of the cart's 3 compartments. The drawers have full-length extension glides and integral pulls. The cart comes in a Soft White powder-coated paint finish.



ECP Cart with Drawers and 2 Center Columns (SRSDCECP)



Materials and Finishes

Base structure: **Powder-coated steel**

Top structure: **304 Stainless Steel**

Corner extrusions: **Clear anodized aluminum**

Corner bumpers: **Thermoplastic polymer (ABS plastic)**

FlexCell panels: **Thermoplastic polymer (ABS plastic)**

Drawers: **Powder-coated steel**

Casters: **Nylon and steel**

Dimensions (nominal) and Weight

Exterior: **57"w x 28.75"d x 36"h**

Interior: **14.5" w (each compartment) x 23.75" d x 12.75"h**

Weight: **330 lbs**

Interior Configuration

- 24 FlexCells, 8 per compartment
- 2 full-width drawers

Options

- Interior and side-mounted exterior accessories



Storage solutions for healthcare

Solairé Medical | P.O. Box 2808 | Grand Rapids, Michigan 49501
888-435-2256 phone | 616-743-5944 fax | solairomedical.com