


Evolve TEE Probe Cabinet with Hinged Door

Evolve modular cabinets create a flexible storage system that organizes a hospital's most supply-intensive departments. A clean design and tight construction offer a refined aesthetic and support infection prevention efforts. The 19.25"-deep, 92"-high TEE probe cabinet comes in 2 widths. It has a solid-

front or tempered glass door, and includes a removable drip tray. Material choices include compression board with a melamine or specified laminate surface, stainless steel, or nonporous AireCore with Glacier White finish, Brushed Aluminum finish, or specified laminate surface. Cabinet has a 4"-high toe kick.



Evolve TEE Probe Cabinet with Hinged Door



Dimensions (nominal)

Exterior (interior):

26.75"w x 19.25"d x 92"h (25.25"w x 17"d x 86.5"h)

36"w x 19.25"d x 92"h (34.5"w x 17"d x 86.5"h) (shown)

92"h cabinet with 6"h toe kick is 94"h

Interior Configuration

26.75"w cabinet holds 3 TEE probes

36"w cabinet holds 5 TEE probes


Options

- Door can be hinged on right or left side of 26.75"w cabinet
- Fan Ventilation Unit
- Keyed, InterConnect lock, and proximity lock upgrade for InterConnect lock
- Sloped top
- 6"h toe kick
- InnerSpace Cloud for remote management and monitoring of cabinet security



Product numbering system for Evolve TEE Probe Cabinet with Hinged Door

Product Shown:
SEC1936PGT



Product Numbers

26.75"w x 19.25"d x 92"h

Compression board: SEC1927PGRT (shown), SEC1927PGLT, SEC1927PSRT, SEC1927PSLT

AireCore: SEA1927PGRT, SEA1927PGLT, SEA1927PSRT, SEA1927PSLT

AireCore and Brushed Aluminum: SEB1927PGRT, SEB1927PGLT, SEB1927PSRT, SEB1927PSLT

Stainless steel: SES1927PGRT, SES1927PGLT, SES1927PSRT, SES1927PSLT



Product Numbers

36"w x 19.25"d x 92"h

Compression board: SEC1936PGT (shown), SEC1936PST

AireCore: SEA1936PGT, SEA1936PST

AireCore and Brushed Aluminum: SEB1936PGT, SEB1936PST

Stainless steel: SES1936PGT, SES1936PST

Evolve TEE Probe Cabinet with Hinged Door

Materials and finishes

Cabinet shell:

- Compression board: Thermally fused melamine over 0.75"-thick compression board or specified Wilsonart, Formica, or Pionite standard laminate
- AireCore: Aluminum composite with Glacier White finish or specified Wilsonart, Formica, or Pionite standard laminate
- AireCore and Brushed Aluminum finish: Aluminum composite with Brushed Aluminum finish
- Stainless steel: 18-gauge 304 stainless steel with 18% chrome and 8% nickel analysis; exposed surfaces

Edge banding:

- Compression board and AireCore: White ABS
- AireCore with Brushed Aluminum finish: Brushed Aluminum acrylic

FlexCell panel:

- Injection-molded ABS

Solid-front door:

- Compression board: Thermally fused melamine over 0.75"-thick compression board or specified Wilsonart, Formica, or Pionite standard laminate
- AireCore: Aluminum composite with Glacier White finish or specified Wilsonart, Formica, or Pionite standard laminate
- AireCore and Brushed Aluminum finish: Aluminum composite with Brushed Aluminum finish
- Stainless steel: 18-gauge 304 stainless steel with 18% chrome and 8% nickel analysis

Glass-front door:

- Compression board: Thermally fused melamine over 0.75"-thick compression board or specified Wilsonart, Formica, or Pionite standard laminate with tempered glass and clear anodized aluminum trim
- AireCore: Aluminum composite with Glacier White finish or specified Wilsonart, Formica, or Pionite standard laminate with tempered glass and clear anodized aluminum trim
- AireCore and Brushed Aluminum finish: Aluminum composite with Brushed Aluminum finish with tempered glass and clear anodized aluminum trim
- Stainless steel: 18-gauge 304 stainless steel with 18% chrome and 8% nickel analysis and tempered glass

Hinges:

- Compression board, AireCore, and AireCore with Brushed Aluminum cabinet: Metal concealed European-style
- Stainless steel cabinet: Stainless steel

Door pulls:

- Compression board, AireCore, and AireCore with Brushed Aluminum cabinet: Dull chrome-plated steel bow pull
- Stainless steel cabinet: 5" stainless steel bow pull